VOCO

Vococid

CE 0482

GB

Etching liquid/gel

Instructions for use

Fields of application:
Vococid is intended to etch only the part of the enamel and/or the dentine (**in total-etch-technique**), which will subsequently be covered by a composite material:
– for composite fillings
– for adhesive fixing of Maryland bridges, inlays, brackets and temporary splinting
– for sealing as caries prophylaxis
– for linings to create a fixed bond (according to manufacturer's instructions)

Indications for the application:
Vococid is both simple to use and easily visible on the enamel due to the blue colouring and consistency of the preparation. With **Vococid** Gel an exact application is possible (e. g. for etching the enamel margins) and the liquid is suitable for treatment of larger surfaces. It is best applied using **Pele Tim** foam pellets or a brush.

Without total-etching (only etching of enamel):
Coat the surface of the enamel, which must be completely moist, with one or two applications of etching liquid/gel on a **Pele Tim**. If necessary, apply extra etching liquid/gel.
Vococid should be allowed to act for the recommended period of time as shown below:
approx. 20 s for normal enamel and for enamel of adolescents;
approx. 30 s for devitalized teeth, deeply fluoridated enamel or enamel without prisms (e.g. primary teeth).
After application rinse thoroughly with water, spraying each etched tooth for at least 15 - 30 s. Afterwards the patient should not close his/her mouth. Dry the etched surfaces thoroughly using an air blower.
Absolutely avoid contamination with saliva. In case of contamination, repeat the etching process (10 s). Check the etched areas (chalky-white).

Total-etch-technique with dentine and enamel bonding materials (etching of enamel and dentine):
Etch entire cavity (enamel appr. 20 - 30 s, dentine max. 15 s). Remove **Vococid** etching gel/liquid by suction, rinse and remove excess moisture with an air-blower.
With other preparations, please refer to the manufacturer's instructions for use.

Warnings/precautions:
Caution! **Vococid** contains phosphoric acid (35%).
In case of contact with the eyes, rinse thoroughly with water. Avoid contact with gingiva, certain composite filling materials and neighbouring teeth.
If total etching is not used, cover the dentine before etching, if necessary. It is recommended to cover adjacent teeth and etched enamel surfaces, which are not covered with composite material, once or twice with a fluoride varnish for remineralisation.
Store at temperatures between 4°C - 23°C. Tightly close after use.
Only give completely emptied packages to recycling centers. **Vococid** residues can be diluted with water or disposed of by neutralising with waste water.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA1063 E1 0310 99 © by VOCO
VOCO GmbH